

# CALIBRATION AND CERTIFICATION SERVICE SOLUTIONS FOR YOUR FACILITY



## THERE WHEN YOU NEED US

EAGLE IS A U.S. FDA AND DEA REGISTERED A2LA ISO 17025 ACCREDITED LABORATORY FOUNDED IN 2003 THAT PROVIDES ANALYTICAL CHEMISTRY AND MICROBIOLOGICAL TESTING, CONSULTATION SERVICES, AND OTHER TECHNICAL SERVICES FOR PHARMACEUTICAL MANUFACTURERS, PHARMACIES, MEDICAL DEVICE COMPANIES, AND OTHER HIGHLY REGULATED INDUSTRIES, TO SUPPORT THEIR COMPLIANCE NEEDS.

OUR ENGINEERS CAN PROVIDE INITIAL AND ONGOING ISO CERTIFICATION OF YOUR FACILITY TO SUPPORT ASEPTIC OPERATIONS, AS WELL AS QUALIFY AND CALIBRATE YOUR EQUIPMENT TO ASSURE PROPER OPERATION AND PERFORMANCE.

---

# OUR SERVICE AND PRODUCT OFFERINGS

**EAGLE HAS EXTENSIVE EXPERIENCE IN CALIBRATION AND CERTIFICATION SOLUTIONS**



Our engineers know that you rely on your equipment to perform as expected. These instruments operate as a safeguard for the quality of your products and to the wellbeing of your employees. To ensure your equipment is operating properly, international standards recommend that most instruments be tested and calibrated.

**With excellent technical and customer services and quick turnaround times, Eagle's engineers will ensure your systems are performing as designed. With Eagle's engineers, you will feel confident that your equipment is operating as intended.**

In addition to calibration and certification services, Eagle offers services in consulting, laboratory testing, validation, and environmental monitoring.



## CERTIFICATION, CALIBRATION AND QUALIFICATION SERVICES

# CERTIFICATION, CALIBRATION, QUALIFICATION

### **CERTIFICATION**

Laminar Airflow Workspaces (LAFW)

Biological Safety Cabinets (BSC)

Powder Containment Hoods

Fume Hoods (coming soon)

Full Cleanroom Certification to <797> & CETA Standards includes:

- HEPA Filter Integrity Testing (Leak Test)

- Differential Pressure Measurement

- Air Change Per Hour Measurement and Calculation

- Non-Viable Air Sampling (Particle Counts)

- Viable Air and Surface Sampling (Environmental Monitoring)

- Static Air Visualization (Smoke) Study

Dynamic Air Visualization (Smoke) Study with Video and Report

### **CALIBRATION**

Balances

Temperature Sensors and Thermometers

Relative Humidity Sensors

Differential Pressure Sensors (digital and magnehelic)

Controlled Temperature Units

Pipettes (coming soon)

### **QUALIFICATION<sup>1</sup>**

Installation/Operational Qualification and Performance Qualification:

- Autoclave - Sterilization Cycle Validation

- Dry-Heat Oven - Depyrogenation Cycle Validation

- Incubators


- Refrigerators

- Freezers

- Stability Chambers

<sup>1</sup> Other equipment available upon request.





# ENVIRONMENTAL MONITORING SOLUTIONS

**TRYPTIC SOY AGAR PLATES - 10 PLATES PER PACK, MINIMUM ORDER - 5 PACKS**

Our tryptic soy agar (TSA) media plates are supplemented with neutralizers to minimize the effect of antimicrobial agents, recovering a wide array of microorganisms.

## **65 MM TSA CONTACT PLATES**

Suitable for active air, surface, and personnel monitoring.

## **90 MM TSA MEDIA PLATES**

Suitable for air and personnel monitoring.

### **FEATURES**

- Room Temperature Storage
- Self-Locking
- Triple Wrapped
- Gamma Irradiated
- Growth Promotion Testing Included
- Quick and Convenient Ordering Process
- Standing Orders for Regular Shipment

### **RELATED SERVICES**

- Incubation and Enumeration of Media
- Microbial Identification
- Temperature Mapping
- Recovery Studies
- Process Control Charts
- Eagle offers on-site training from our microbiologists to your personnel so you can have documented training on the evaluation of plates after incubation.

For additional information, visit [eagleanalytical.com/em](http://eagleanalytical.com/em).





## ENVIRONMENTAL MONITORING SOLUTIONS

### SURFACESHIELD WIPE SAMPLING AND CLEANING VERIFICATION KIT

The SurfaceShield Wipe Sampling and Cleaning Verification Kit is a fast, easy, and convenient method for verifying the efficacy of your cleaning program, including deactivation, decontamination, cleaning, and/or disinfection for sterile. This multipurpose and affordable kit can be utilized to achieve and maintain USP <800> compliance and meet Title 21 Part 211 CFR requirements.

#### BASE KIT INCLUDES

- 2 - Sample Sets (tube, swab, single-use vial of sterile water, pair of gloves in reclosable plastic bag, 10 x 10 cm (100 cm<sup>2</sup>) template)
- 3 oz single-use cold pack
- Sampling Instructions insert
- Blank labels
- Insulated mailer

For pricing visit [eagleanalytical.com/testing-services-pricing](http://eagleanalytical.com/testing-services-pricing).



---

# CONTACT US

WE ARE HERE TO HELP



Visit our [website](https://www.eagleanalytical.com) or scan the QR code to learn more about our team, calibration and certification capabilities, or to receive an initial consultation to discuss your operation's specific needs. We look forward to assisting you.



[WWW.EAGLEANALYTICAL.COM](https://www.eagleanalytical.com)

1.800.745.8916 | [INFO@EAGLEANALYTICAL.COM](mailto:INFO@EAGLEANALYTICAL.COM) | 11111 S. WILCREST DR. S1000 HOUSTON, TEXAS 77099

